



STOPPERSOCKEN

REF 945030

Product Description, Intended use, Application

The **size XXL** grey sock is non-slip thanks to the patterned anti-slip pads on the top and tread of the sock. When worn, the product offers the patient more sure-footedness and traction when walking, so that falls and subsequent injuries can be reduced. Stopper studs on both sides ensure safe walking even if the sock are twisted. The socks are individually packed in pairs.

Composition

Cotton, polyester elastane, polyvinylchloride

Contra-Indications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: film bag

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.